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
Organ / System	Disease model / Process	Target	Formulation / Delivery Reagent	Injection / Delivery	Dosing / kg animal / injection	Volume/injection	Frequency	Reference
Brain	Alzheimer	APP	None	Intracerebral	8.5 mg	6 µl/day	4 weeks	Poster on IVGN website
	ADHD	Dopamine	None	Intracerebral	20 mg	6 µl to 12 µl/day	1 or 2 weeks	Thakker et al. 2004
		5-HT transporter	None	Intracerebral	20 mg	6 µl to 12 µl/day	1 or 2 weeks	Thakker et al. 2005
	Depression	Serotonin transporter	None	Intracerebral	20 mg	6 µl to 12 µl/day	1 or 2 weeks	Thakker et al. 2005
	Glutamate release	GluR2, COx-1	None	Intracerebral	N/A	N/A	Once	Akaneya et al. 2005
Eye	Ocular Inflammation and fibrosis	TGFβ RII	Liposome	Local	0.005	30 µl	Once	Nakamura et al. 2004
	Herpetic Stromal keratitis	VEGF/R	None	Local	0.5 mg	10 µl	at 6 hrs and 24 hrs	Kim et al. 2004
		VEGF/R	Pegylated	Intravenous	2 mg	100 µl	at 6 hrs and	Kim et al. 2004

		R	immunoliposomes				24 hrs	
Hearing	Hearing-autosomal dominant	GJB2	Liposome	Local	12.5 µg-50 µg	1-2 µl	Applied to gel foam, placed locally	Maeda et al. 2005
	Rheumatoid arthritis	TNFα	None	Local	0.5 mg	N/A	weekly (for 3 weeks)	Schiffelers et al. 2005
Inflammation		TNFα	PEI	Intravenous	2 mg	N/A	Once	Schiffelers et al. 2005
		IL-12p40	Oligofectamine™	Intraperitoneal	N/A	200 µl	N/A	Flynn et al. 2004
Spinal cord	Neuropathic pain	P2X3 cation channel	None	Intrathecal	20 mg	6 to 12 µl/day	Every day for 6-7 days	Dorn et al. 2004
	Phrenic long term facilitation	BDNF	Oligofectamine™	Intracerebral	0.6 mg	20 µl	Once	Baker-Herman et al. 2004
	Neuropathic pain	DOR	Cationic lipid	Intrathecal	0.1 mg	10 µl	Daily at 8 hrs and	Luo et al. 2005
	Acute liver Failure	Fas	None	Hydrodynamic	2.5 mg	1 ml	24 hrs	Song et al. 2003
Liver		Caspase 8	None	Hydrodynamic	0.325 mg	2.5 ml (tail Vein)	Once	Zender et al. 2003
		Caspase 8	None	Hydrodynamic	0.325 mg	1 ml portal vein	Once	Zender et al. 2003
	Ischemia/reperfusion	Caspase 8/3	10% lipiodol	Hydrodynamic	0.375 mg	1 ml portal Vein	Once	Contreras et al. 2004
	Hypercholesterolemia	ApoB	Modified-Cholesterol conjugate	Intravenous	1 mg	200 µl	Over 3 days	Soutschek et al. 2004
Kidney	Renal ischemia/reperfusion	Fas	None	Hydrodynamic	2.5 mg	1 ml	N/A	hamar et al. 2004
Lung	Ischemia/reperfusion	Heme oxygenase1	None	Intranasal	2 mg	50 µl	Once	Zhang et al. 2004
Metabolism	Obesity	AGRP	None	Intracerebral	0.2 mg	0.5 µl	Once	Makimura et al. 2002
Pancrea	Diabetes	Ins2	Lipofectamine™	Local injection	0.2 mg	100 µl	N/A	Bradley et al. 2005

s		2000						
		Ins2	PBS	Intravenous	5 mg	800 µl	N/A	Bradley et al. 2005
Bone cancer		EZH2, p110-a	atelocollagen	Intravenous	2.5 mg	200 µl	At 3, 6 and 9 Days	Takeshita et al. 2005
Breast cancer		CXCR 4	None	Intravenous	0.5 mg	N/A	Twice weekly	Liang et al. 2005
		Colony - stimulating factor	None	Local	0.5 mg	N/A	Every 3 days	Aharinejad et al. 2004
		RhoA	None	Local	6.75 µg	100 µl	Every 3 days (for 20 days)	Pille et al. 2005
		c-raf	liposome	Intravenous	7.5 mg	N/A	2 injections a day for 5 days	Chein et al.
		Raf-1	Polymere	Local	0.2 mg	N/A	Every 5 days	Leng et al. 2005
Cancer	Ewing's sarcoma	VEGF	PEI	Local	20 µg	N/A	Twice weekly	Guan et al. 2005
	Germ cell tumor	FGF4	atelocollagen	Local	0.125 mg	50 µl	Once	Minakuchi et al. 2004
	Liver cancer	bcl-2	Liposome	Intravenous	10 mg	N/A	Two 5-day cycles with daily injection	Yano et al. 2004
		Spingosine-1 phosphate-R	Liposome	Local	0.75 mg	60 µl	Every 3 days	Chae et al. 2004
	Lung cancer / Cervical cancer	Polo-like kinase 1	None (ATA-treated)	Intravenous	0.5 mg	500 µl PBS	3 times per week	Spankuch et al. 2004
	Melanoma	B-RAF	Ex vivo, electroporation	subcutaneous	N/A	N/A	NA	Sharma et al. 2005

	ErbB2	Oligofectamine	Intratumoral, Intravenous	2 mg	N/A	N/A	Song et al. 2005
Neuroblastoma	VEGF-R erbb2/n eu (HER2) )R	Ligand targeted	Intravenous	2 mg	200 µl	Every 3 Days	Schiffelers et al. 2004
Ovarian cancer		PEI-complex	Intraperitoneal	0.45 mg	N/A	Every 2-3 days	Urban-klein et al. 2005
Pancreatic cancer	bcl-2	None	Intraperitoneal	0.2 mg	120 µl	Daily for 24 Days	Ocker et al. 2005
	FAK	None	Intravenous	0.15 mg	50 µl	Twice weekly	Duxbury et al. 2003
	Epha2	None	Intravenous	0.15 mg	50 µl	Twice weekly	Duxbury et al. 2004
	CEAC AM6	None	Intravenous	0.15 mg	N/A	Twice weekly	Duxbury et al. 2004
Prostate cancer	EphB4	None	Intraperitoneal	10 mg	N/A	Daily for 2 weeks	Xia et al.
	VEGF	atelocollagen	Local./Intravenous	0.35 mg	50 µl	Every 10 Days	Takei et al. 2004
West Nile virus	viral genes	None	Hydrodynamic	9 mg	1.6 ml Ringer sol	Once	Bai et al. 2005
Hepatitis C virus	viral genes	None (use RNasin)	Hydrodynamic	2 mg	1.8 ml PBS	N/A	McCaffrey et al. 2002
Coxsackievirus B3	viral genes	None	Hydrodynamic	9 mg	1.6 ml NACl 0.9%	Once	Merl et al. 2005
RSV, Parainfluenza	viral genes	None/Liposome	Intranasal	3.5 mg	41 µl	N/A	Bitko et al. 2005
Influenza	viral genes	PEI	Intravenous/ Intranasal	3 mg	200 µl	Once	Ge et al. 2004
Infectious agents	viral genes	None	Hydrodynamic	2.8 mg	1 ml	Once	Tompkins et al. 2004
	viral genes	Oligofectamine™	Intranasal	1.15 mg	50 µl	Once	Tompkins et al. 2004
Hepatitis B	HBsAg	None	Hydrodynamic	1.25 mg	1.5 ml	Once	Giladi H et al. 2003
	HBsAg	None	Hydrodynamic	N/A	N/A	N/A	Klein et al. 2003

	HBsAg	None	Hydrodynamic	0.5 mg	2 ml	Once	McCaffrey et al. 2003
	Viral genes	SNALP	Intravenous	1,3 or 5 mg	N/A	N/A	Morrissey et al. 2005
Herpes simplex virus 2 (HSV-2)	Viral genes	Oligofectamine™	Intravaginal	0.36 mg	12 µl	Twice	Palliser et al. 2005
Sepsis	TNFα	None	Intraperitoneal	1.5 mg	N/A	Once	Sorensen et al. 2003

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